



www.CloudNeuro.com



To Get Started Call 877-835-3364
or visit us online at www.CloudNeuro.com

27 Arlington Road, Woburn, MA 01801, USA

STATE OF THE ART RTMS SYSTEM

Start your TMS Practice Today

CloudTMS

A Leader in rTMS since 1996

- CloudTMS has been in the TMS business for over 20 years.
- Our first stimulator was developed in 1996
- Thousands of CloudTMS devices are currently being used in Neurological and Psychiatric clinics around the World

Main Unit

- Controls System operations
- Can also be used in Stand-Alone mode
- Full Front Panel Control and Display of all Parameters including Coil Temperature
- USB Connection to Laptop

Intended Use

The CloudTMS device is indicated for the treatment of Major Depressive Disorders in adult patients who have failed to receive satisfactory improvement from prior antidepressant medications in the current episode.

Smart Windows Based Software: Touchscreen Optimized

- Windows-based software for Preset Treatment Protocols
- Easy to configure to fit your Practice/Patient needs
- Sets follow-up Sessions Stimulator Parameters and Coil Positioning hints
- Optimized for touchscreen interface

Special Cooling Unit: Cooled Coils

- Cooling liquid runs through the Coil and the Cooler
- Most effective Synthetic Oil cooling system in the Industry
- Makes Overheating a thing of the past





Extra Power Supply Unit

- Main unit designed to operate at up to 30 Hz
- Maximum intensity reached at 5-7 Hz
- Extra Power Supply Unit makes it possible to sustain stimulation Frequencies up to 20-25 Hz, and can deliver a maximum Frequency of up to 100 Hz

Coil Holder

Your special coil holder makes it a breeze to freely move the coil and quickly and safely position it.

Coils

 1.1 T Maximal induction	Cooled	Non-Cooled
	Figure-of-eight coil FEC-02-100-C	Figure-of-eight coil FEC-02-100
 1.6 T Maximal induction	Cooled	Non-Cooled
	Figure-of-eight coil AFEC-02-100-C	Figure-of-eight coil AFEC-02-100